

Foniatr - Feature #240

U stimulu umožnit definovat interval, který se ma prehrat

02 July 2015 13:19 - Tomáš Kulhánek

Status:	Resolved	Start date:	25 June 2015
Priority:	Normal	Due date:	
Assignee:	Tomáš Kulhánek	% Done:	0%
Category:	HTML5 client	Estimated time:	8.00 hours
Target version:	Psychotest2	Spent time:	8.00 hours
Description			
napr 1.wav(4.1-6.4s) tj od 4.1 do 6.4 sekundy			
Subtasks:			
Task # 233: Novy typ tasku - hodnotit segmenty zvuku			New

History

#1 - 09 July 2015 10:07 - Tomáš Kulhánek

- Status changed from New to Assigned

pridat nastroj pro prehrani stimulu a urceni intervalu

#2 - 10 July 2015 11:37 - Tomáš Kulhánek

- Category set to HTML5 client

- Status changed from Assigned to Resolved

- % Done changed from 0 to 100

Implementovano

Interval se definuje za nazvem zvuku v hranatých závorkách ve tvaru [starttime,duration] v milisekundách

stimulus trombone.wav [1000,3000]

prehraje od 1 sekundy po dobu 3 sekundy. Test viz napr.

<http://physiome.lf1.cuni.cz/psychoacoustictest/psychotest.htm#119>